



PODŁOKIETNIK

PRZEZNACZENIE

1. Dla osób ze spastyką oraz przykurczami dłoni i nadgarstka.


DZIAŁANIE

1. Dopasowanie do anatomicznego kształtu ręki.
2. Odwiedzenie i wyprostowanie palców i nadgarstka.
3. Stabilne podparcie ramienia ma bezpośredni wpływ na prawidłowe ustawienie kręgosłupa.

ZASTOSOWANIE

1. Prawidłowe ułożenie i ustabilizowanie przedramienia przy hemiplegii.

WYMIARY



ZAGŁÓWEK

PRZEZNACZENIE

1. Stabilizacja dla osób z obniżoną kontrolą głowy.


DZIAŁANIE

1. Dopasowuje się do kształtu głowy.
2. Umożliwia podparcie głowy w najbardziej korzystnym punkcie.

ZASTOSOWANIE

1. Uniwersalny zagłówek do większości wózków.

WYMIARY



VELCROARMREST



INDICATIONS

1. For persons with spasticity and hand and wrist contractures.

ACTION

1. Adjusted to the anatomical shape of the hand.
2. Separating and straightening of the fingers and wrist.
3. A stable support for the arm has a direct influence on the proper positioning of the spine.

USE

1. Proper positioning and stabilization of the forearm in the case of hemiplegia.

MEASUREMENTS

[cm]	A	B
S	25	30
M	25	45

HEADREST



INDICATIONS

1. Stabilization for persons with reduced head control.

ACTION

1. Adapts to the shape of the head.
2. Provides support for the head at the most-convenient point.

USE

1. A universal headrest for most wheelchairs.

MEASUREMENTS

[cm]	A	B
M	15	30

mount 32 or 42
mocowanie 32 lub 42